



Oxyrase, Inc.
3000 Park Ave. West
Mansfield, OH 44906
Ph.: 419-589-8800
Fax: 419-589-9919
www.oxyrase.com

Oxyrase CryoBroth™ Product Insert

Oxyrase CryoBroth™ is a sterile enzyme formulation that provides a protective anaerobic environment by utilizing Nature's Antioxidant® for preservation and enhanced recovery of organisms that have been stored cryogenically.

Precautions:

Oxyrase CryoBroth™ is a filter sterilized product and must be handled aseptically to maintain sterility. A **Safety Data Sheet** is available on our website.

Product Characteristics:

The formulation of Oxyrase CryoBroth™ combines cryo-protectants (to minimize damage from freezing / thawing) and Nature's Antioxidant® (to enhance the recovery of cells from the frozen state).

Freezing and thawing causes cell damage (3), Oxyrase CryoBroth™ has been used for the storage and recovery of anaerobes (including facultative anaerobes) and numerous other cell cultures (1). Oxyrase® Enzyme System enhances recovery of injured cells (2), and has been used to preserve cells frozen and stored in cryogenic conditions (3,4).

Limitations:

Oxyrase CryoBroth™ contains a penicillin binding protein that may interfere with penicillin and some related antibiotics.

Handling and Storage Instructions:

Oxyrase CryoBroth™ will arrive thawed but cold. Storage options are listed below:

- Long Term Storage:** Store the product at $\leq -20^{\circ}\text{C}$ to maintain full activity. Oxyrase CryoBroth™ can be thawed and re-frozen five times without affecting its activity and performance. In cases where the product will be used infrequently and / or in small amounts, aseptically aliquot the product into smaller, individual, sterile containers (refer to short term storage, if needed).
- Short Term Storage:** Store the product at 2°C to 8°C for use within 30 days (a precipitant may form at this temperature).

When stored in this manner, the product will maintain its full activity to the printed expiration date on the label.

Thawing Oxyrase CryoBroth™:

A convenient way to thaw Oxyrase CryoBroth™ is to place it in the refrigerator overnight.

If necessary, the product can be thawed by warming. *Do not*

exceed a warming temperature of 37°C . Only apply heat to the outside of the container while ice is still present inside the container. When all ice has melted, keep the product chilled by placing the container in ice until ready for use.

To ensure uniform activity within a thawed sample, *gently* mix the product before use or distribution (*do not agitate vigorously*). Vigorous agitation (i.e. shaking) causes foaming and denatures protein in the product, which may result in loss of activity.

In some cases, precipitate may be observed, but will not affect Oxyrase CryoBroth™ performance.

Instructions for Use:

Aseptically, aliquot Oxyrase CryoBroth™ volume into a sterile container. Add a high density suspension of cells ($>10^9$ cfu/mL is recommended for maximum recovery) directly into the Oxyrase CryoBroth™. *Gently* mix (*do not agitate vigorously or aerate the cell suspension*) the suspension using a sterile disposable pipette. Aseptically, add the cell suspension to each container. Place each container in a cryogenic freezer immediately after distribution.

Alternately, make, mix, and distribute 0.2 mL of the above cell suspension into cryo-vials.

When thawing, allow 15 - 30 minutes at room temperature for vial(s) to thaw completely. Immediately, transfer cells to the growth medium appropriate for organism cultivation.

As a control, it is recommended that a single vial be thawed 7-14 days after freezing to ensure / test viability of the cell suspension made. This provides a base-line for measurement of cell recovery, and serves as a point of comparison to measure performance each time the Oxyrase CryoBroth™ cell suspension is reconstituted.

Quality Control:

The length of time that cells can remain viable in Oxyrase CryoBroth™ is dependent on numerous factors: nature, type of cell culture, and storage conditions of cryo-vials.

It is recommended that recovery of Oxyrase CryoBroth™ cell suspensions be scored and recorded each time that an Oxyrase CryoBroth™ vial is used. This history of cell viability, along with storage time, provides a useful record to determine optimal storage conditions.

Guarantee:

Oxyrase CryoBroth™ has a shelf-life of 12 months from date of manufacture, under recommended storage and use conditions.

If Oxyrase CryoBroth™ does not preserve and enhance recovery of organisms that have been stored cryogenically as specified under recommended storage and use conditions, Oxyrase, Inc. will refund your purchase price. To receive a product refund, write or call Oxyrase Inc. with the product lot number which is located on the Oxyrase CryoBroth™ label. Oxyrase, Inc. is available to answer any questions about this product and its applications.

Oxyrase, Inc. LAB.0023.v.004 CRF#172

- Linda S.L. Yu and Daniel Y., C. Fung. 1990. Effect of Oxyrase Enzyme on *Listeria Monocytogenes* and Other Facilitative Anaerobes. *J. Food Safety* 11:3, 163-175.
- SM George, LLC Richardson, IE Pol, MW Peck. 1998. Effect of Oxygen Concentration and Redox potential on Recovery of Sub-Lethally Heat-Damaged Cells of *E. coli* 0157:H7, *Salmonella enteritidis* and *Listeria monocytogenes*. *Appl. Microbiol.* 84:5, 903-909.
- Mazur, Kathov, Katkova, and Critser. 2000. The Enhancement of the Ability of Mouse Sperm to Survive Freezing and Thawing by the Use of High Concentrations of Glycerol and Presence of an *Escherichia coli* Membrane Preparation (Oxyrase) to Lower the Oxygen Concentration. *Cryobiology.* 40:4, 187-209.
- Koshimoto, Gamliel, Mazur. 2000. Effect of Osmolality and Oxygen Tension on the Survival of Mouse Sperm Frozen to Various Concentrations of Glycerol and Raffinose. *Cryobiology.* 41:3, 204-231.