


Oxidation Test *Porphyromonas levii*



Non-PRAS

Post Reduced

PRAS Brucella

PRAS 2

PRAS 3


PRAS Schaedler

Oxidation of Post-Reduced Blood Agar Plates (BAP)

Oxyrase PRAS Enzymatic Reduction **10:00am Start** Competitor Non-PRAS Post-Reduced



10:40am



24 hours



In terms of color change as a sign of oxidation, the enzymatically reduced PRAS plate began to oxidize in approximately 24 hours. Standard Post-Reduced BAP plate oxidized in 40 minutes.